



LENS AVAILABILITY + RANGES

STOCK LENS OPTIONS + RANGES



14505 N Hayden Rd
Suite #101
Scottsdale AZ 85260
(800) 833-4779

© May 2024 GSRx®, Inc.

www.GSRx.com

ND⁴™ & epik™ Variable Progressive • Ultra High Definition Vision
Use for patients demanding the best all around performance.

Index	Material Name	ADD	Power Ranges	CYL
1.50	CR39 Plastic	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.50	CR39 Transitions®	ADD .75 to 3.50	-9.50 to +7.50	up to -6
1.50	CR39 Polarized	ADD .75 to 3.50	-9.25 to +7.00	up to -6
1.50	CR39 XtrActive	ADD .75 to 3.50	-11.75 to +6.25	up to -6
1.50	CR39 epik™ BLUE	ADD .75 to 3.50	-10.50 to +8.00	up to -6
1.53	Trivex®	ADD .75 to 3.50	-10.00 to +7.00	up to -6
1.53	Trivex® Transitions®	ADD .75 to 3.50	-7.50 to +5.50	up to -4
1.53	Trivex® Polarized	ADD .75 to 3.50	-7.00 to +7.00	up to -4
1.53	Trivex® epik™ BLUE	ADD .75 to 3.50	-8.00 to +6.00	up to -5
1.59	PolyTrue™	ADD .75 to 3.50	-11.00 to +7.00	up to -6
1.59	Poly Transitions®	ADD .75 to 3.50	-10.50 to +8.00	up to -6
1.59	Poly Polarized	ADD .75 to 3.50	-10.75 to +7.00	up to -6
1.59	Poly epik™ BLUE	ADD .75 to 3.50	-11.00 to +7.00	up to -6
1.60	1.60 High Index	ADD .75 to 3.50	-12.50 to +9.50	up to -6
1.60	1.60 Transitions®	ADD .75 to 3.50	-12.00 to +8.50	up to -6
1.60	1.60 Polarized	ADD .75 to 3.50	-11.50 to +7.00	up to -6
1.60	1.60 XtrActive	ADD .75 to 3.50	-11.00 to +7.50	up to -6
1.60	1.60 epik™ BLUE	ADD .75 to 3.50	-11.00 to +9.00	up to -6
1.67	1.67 High Index	ADD .75 to 3.50	-13.00 to +12.50	up to -6
1.67	1.67 Transitions®	ADD .75 to 3.50	-11.50 to +8.00	up to -6
1.67	1.67 Polarized	ADD .75 to 3.50	-12.00 to +8.50	up to -6
1.67	1.67 XtrActive	ADD .75 to 3.50	-12.00 to +8.50	up to -6
1.67	1.67 epik™ BLUE	ADD .75 to 3.50	-12.50 to +11.00	up to -6
1.74	1.74 High Index	ADD .75 to 3.50	-18.00 to +7.50	up to -6
1.74	1.74 Transitions®	ADD .75 to 3.50	-17.50 to +11.00	up to -6
1.74	1.74 epik™ BLUE	ADD .75 to 3.50	-18.00 to +7.50	up to -6

Fitting heights variable 13-35mm; fitting cross 4mm.

This lens design is optimized for patient Rx, is Rx Compensated, and has POW capabilities. Power ranges vary due to blank thickness and manufacturer. May also be affected by patient frame and Rx (prism).

ND⁴™ & epik™ Variable Device Progressive • Ultra High Definition Vision
Use for patients who are high device users (close in).

Index	Material Name	ADD	Power Ranges	CYL
1.50	CR39 Plastic	ADD .75 to 3.50	-9.25 to +7.50	up to -6
1.50	CR39 Transitions®	ADD .75 to 3.50	-9.50 to +7.50	up to -6
1.50	CR39 Polarized	ADD .75 to 3.50	-6.50 to +5.50	up to -6
1.50	CR39 epik™ BLUE	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.53	Trivex®	ADD .75 to 3.50	-9.00 to +6.00	up to -6
1.53	Trivex® Transitions®	ADD .75 to 3.50	-7.50 to +5.00	up to -4
1.53	Trivex™ Polarized	ADD .75 to 3.50	-4.50 to +6.00	up to -4
1.53	Trivex® epik™ BLUE	ADD .75 to 3.50	-8.00 to +6.00	up to -5
1.59	PolyTrue™	ADD .75 to 3.50	-9.00 to +7.00	up to -6
1.59	Poly Transitions®	ADD .75 to 3.50	-9.50 to +7.50	up to -6
1.59	Poly Polarized	ADD .75 to 3.50	-10.00 to +7.00	up to -6
1.59	Poly epik™ BLUE	ADD .75 to 3.50	-9.00 to +7.00	up to -6
1.60	1.60 High Index	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.60	1.60 Transitions®	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.60	1.60 Polarized	ADD .75 to 3.50	-10.00 to +7.00	up to -6
1.60	1.60 epik™ BLUE	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.67	1.67 High Index	ADD .75 to 3.50	-12.00 to +8.00	up to -6
1.67	1.67 Transitions®	ADD .75 to 3.50	-12.00 to +8.00	up to -6
1.67	1.67 Polarized	ADD .75 to 3.50	-12.00 to +8.00	up to -6
1.67	1.67 epik™ BLUE	ADD .75 to 3.50	-11.50 to +8.00	up to -6
1.74	1.74 High Index	ADD .75 to 3.50	-15.00 to +7.50	up to -6
1.74	1.74 Transitions®	ADD .75 to 3.50	-15.00 to +7.50	up to -6
1.74	1.74 epik™ BLUE	ADD .75 to 3.50	-15.00 to +7.50	up to -6

Fitting heights variable 13-35mm; fitting cross 4mm.

This lens design is optimized for patient Rx, is Rx Compensated, and has POW capabilities. Power ranges vary due to blank thickness and manufacturer. May also be affected by patient frame and Rx (prism).

ND4™ & epik™ Variable Active Progressive

Ultra High Definition Vision • Use for patients who participate in active sports.

Index	Material Name	ADD	Power Ranges	CYL
1.50	CR39 Plastic	ADD .75 to 3.50	-9.25 to +7.50	up to -6
1.50	CR39 Transitions®	ADD .75 to 3.50	-9.50 to +7.50	up to -6
1.50	CR39 Polarized	ADD .75 to 3.50	-6.50 to +5.50	up to -6
1.50	CR39 epik™ BLUE	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.53	Trivex®	ADD .75 to 3.50	-9.00 to +6.00	up to -6
1.53	Trivex® Transitions®	ADD .74 to 3.50	-7.50 to +5.00	up to -4
1.53	Trivex® Polarized	ADD .75 to 3.50	-4.50 to +6.00	up to -4
1.53	Trivex® epik™ BLUE	ADD .75 to 3.50	-8.00 to +6.00	up to -5
1.59	PolyTrue™	ADD .75 to 3.50	-9.00 to +7.00	up to -6
1.59	Poly Transitions®	ADD .75 to 3.50	-9.50 to +7.50	up to -6
1.59	Poly Polarized	ADD .75 to 3.50	-10.00 to +7.00	up to -6
1.59	Poly epik™ BLUE	ADD .75 to 3.50	-9.00 to +7.00	up to -6
1.60	1.60 High Index	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.60	1.60 Transitions®	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.60	1.60 Polarized	ADD .75 to 3.50	-10.00 to +7.00	up to -6
1.60	1.60 epik™ BLUE	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.67	1.67 High Index	ADD .75 to 3.50	-12.00 to +8.00	up to -6
1.67	1.67 Transitions®	ADD .75 to 3.50	-12.00 to +8.00	up to -6
1.67	1.67 Polarized	ADD .75 to 3.50	-12.00 to +8.00	up to -6
1.67	1.67 epik™ BLUE	ADD .75 to 3.50	-11.50 to +8.00	up to -6
1.74	1.74 High Index	ADD .75 to 3.50	-15.00 to +7.50	up to -6
1.74	1.74 Transitions®	ADD .75 to 3.50	-15.00 to +7.50	up to -6
1.74	1.74 epik™ BLUE	ADD .75 to 3.50	-15.00 to +7.50	up to -6

*Fitting heights variable 13-35mm; fitting cross 4mm.**This lens design is optimized for patient Rx, is Rx Compensated, and has POW capabilities. Power ranges vary due to blank thickness and manufacturer. May also be affected by patient frame and Rx (prism).***ND4™ Fixed Progressive • Fixed Corridor Ultra High Definition Vision**

Use for patients demanding the best all around performance.

Index	Material Name	ADD	Power Ranges	CYL
1.50	CR39 Plastic	ADD .75 to 3.50	-9.25 to +7.50	up to -6
1.50	CR39 Transitions®	ADD .75 to 3.50	-9.50 to +7.50	up to -6
1.50	CR39 Polarized	ADD .75 to 3.50	-6.50 to +5.50	up to -6
1.50	CR39 epik™ BLUE	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.53	Trivex®	ADD .75 to 3.50	-9.00 to +6.00	up to -6
1.53	Trivex® Transitions®	ADD .75 to 3.50	-7.50 to +5.00	up to -4
1.53	Trivex® Polarized	ADD .75 to 3.50	-4.50 to +6.00	up to -4
1.53	Trivex® epik™ BLUE	ADD .75 to 3.50	-8.00 to +6.00	up to -5
1.59	PolyTrue™	ADD .75 to 3.50	-9.00 to +7.00	up to -6
1.59	Poly Transitions®	ADD .75 to 3.50	-9.50 to +7.50	up to -6
1.59	Poly Polarized	ADD .75 to 3.50	-10.00 to +7.00	up to -6
1.59	Poly epik™ BLUE	ADD .75 to 3.50	-9.00 to +7.00	up to -6
1.60	1.60 High Index	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.60	1.60 Transitions®	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.60	1.60 Polarized	ADD .75 to 3.50	-10.00 to +7.00	up to -6
1.60	1.60 epik™ BLUE	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.67	1.67 High Index	ADD .75 to 3.50	-12.00 to +8.00	up to -6
1.67	1.67 Transitions®	ADD .75 to 3.50	-12.00 to +8.00	up to -6
1.67	1.67 Polarized	ADD .75 to 3.50	-12.00 to +8.00	up to -6
1.67	1.67 epik™ BLUE	ADD .75 to 3.50	-11.50 to +8.00	up to -6
1.74	1.74 High Index	ADD .75 to 3.50	-15.00 to +7.50	up to -6
1.74	1.74 Transitions®	ADD .75 to 3.50	-15.00 to +7.50	up to -6
1.74	1.74 epik™ BLUE	ADD .75 to 3.50	-15.00 to +7.50	up to -6

*Fitting heights 14mm, 17mm, 20mm; fitting cross 4mm.**This lens design is optimized for patient Rx, is Rx Compensated, and has POW capabilities. Power ranges vary due to blank thickness and manufacturer. May also be affected by patient frame and Rx (prism).*

epik™ Day Progressive

Ultra High Definition Vision • Use for all patients as your premium daily lens.

Index	Material Name	ADD	Power Ranges	CYL
1.50	CR39 Plastic	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.50	CR39 Transitions®	ADD .75 to 3.50	-9.00 to +5.50	up to -6
1.50	CR39 Polarized	ADD .75 to 3.50	-8.00 to +7.00	up to -6
1.50	CR39 epik™ BLUE	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.53	Trivex®	ADD .75 to 3.50	-9.00 to +6.00	up to -6
1.53	Trivex® Transitions®	ADD .75 to 3.50	-7.50 to +5.00	up to -4
1.53	Trivex® Polarized	ADD .75 to 3.50	-4.50 to +6.00	up to -4
1.53	Trivex® epik™ BLUE	ADD .75 to 3.50	-8.00 to +6.00	up to -5
1.59	PolyTrue™	ADD .75 to 3.50	-9.00 to +7.00	up to -6
1.59	Poly Transitions®	ADD .75 to 3.50	-9.50 to +7.50	up to -6
1.59	Poly Polarized	ADD .75 to 3.50	-10.00 to +7.00	up to -6
1.59	Poly epik™ BLUE	ADD .75 to 3.50	-9.00 to +7.00	up to -6
1.60	1.60 High Index	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.60	1.60 Transitions®	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.60	1.60 Polarized	ADD .75 to 3.50	-10.00 to +7.00	up to -6
1.60	1.60 epik™ BLUE	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.67	1.67 High Index	ADD .75 to 3.50	-12.00 to +8.00	up to -6
1.67	1.67 Transitions®	ADD .75 to 3.50	-12.00 to +8.00	up to -6
1.67	1.67 Polarized	ADD .75 to 3.50	-12.00 to +8.00	up to -6
1.67	1.67 epik™ BLUE	ADD .75 to 3.50	-11.50 to +8.00	up to -6
1.74	1.74 High Index	ADD .75 to 3.50	-15.00 to +7.50	up to -6
1.74	1.74 Transitions®	ADD .75 to 3.50	-15.00 to +7.50	up to -6
1.74	1.74 epik™ BLUE	ADD .75 to 3.50	-15.00 to +7.50	up to -6

*Fitting heights 13, 15, 17, 19, 21mm; fitting cross 4mm.**This lens design is optimized for patient Rx, is Rx Compensated, and has POW capabilities. Power ranges vary due to blank thickness and manufacturer. May also be affected by patient frame and Rx (prism).***epik™ Standard Progressive**

Enhanced Definition Vision • Use as an upgrade to regular progressives.

Index	Material Name	ADD	Power Ranges	CYL
1.50	CR39 Plastic	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.50	CR39 Transitions®	ADD .75 to 3.50	-9.00 to +5.50	up to -6
1.50	CR39 Polarized	ADD .75 to 3.50	-8.00 to +7.00	up to -6
1.50	CR39 epik™ BLUE	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.53	Trivex®	ADD .75 to 3.50	-9.00 to +6.00	up to -6
1.53	Trivex® Transitions®	ADD .75 to 3.50	-7.50 to +5.00	up to -4
1.53	Trivex® Polarized	ADD .75 to 3.50	-4.50 to +6.00	up to -4
1.53	Trivex® epik™ BLUE	ADD .75 to 3.50	-8.00 to +6.00	up to -5
1.59	PolyTrue™	ADD .75 to 3.50	-10.00 to +5.50	up to -6
1.59	Poly Transitions®	ADD .75 to 3.50	-10.50 to +8.00	up to -6
1.59	Poly Polarized	ADD .75 to 3.50	-10.00 to +7.00	up to -6
1.59	Poly epik™ BLUE	ADD .75 to 3.50	-11.00 to +8.00	up to -6
1.60	1.60 High Index	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.60	1.60 Transitions®	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.60	1.60 Polarized	ADD .75 to 3.50	-10.00 to +7.00	up to -6
1.60	1.60 epik™ BLUE	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.67	1.67 High Index	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.67	1.67 Transitions®	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.67	1.67 Polarized	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.67	1.67 epik™ BLUE	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.74	1.74 High Index	ADD .75 to 3.50	-10.00 to +7.50	up to -6
1.74	1.74 Transitions®	ADD .75 to 3.50	-10.00 to +7.50	up to -6
1.74	1.74 epik™ BLUE	ADD .75 to 3.50	-10.00 to +7.50	up to -6

*Fitting heights 14mm and 17mm; fitting cross 4mm.**This lens design is optimized for patient Rx, is Rx Compensated, and has POW capabilities. Power ranges vary due to blank thickness and manufacturer. May also be affected by patient frame and Rx (prism).*

epik™ Office UltraHD

Book, Desk, Room

Index	Material Name	ADD	Power Ranges	CYL
1.50	CR39 Plastic	ADD 1.50 to 3.50	-9.25 to +7.50	up to -6
1.50	CR39 Transitions®	ADD 1.50 to 3.50	-9.50 to +7.50	up to -6
1.50	CR39 epik™ BLUE	ADD 1.50 to 3.50	-10.00 to +8.00	up to -6
1.53	Trivex®	ADD 1.50 to 3.50	-9.00 to +6.00	up to -6
1.53	Trivex® Transitions®	ADD 1.50 to 3.50	-7.50 to +5.00	up to -4
1.53	Trivex® epik™ BLUE	ADD 1.50 to 3.50	-8.00 to +6.00	up to -5
1.59	PolyTrue™	ADD 1.50 to 3.50	-9.00 to +7.00	up to -6
1.59	Poly Transitions®	ADD 1.50 to 3.50	-9.50 to +7.50	up to -6
1.59	Poly epik™ BLUE	ADD 1.50 to 3.50	-9.00 to +7.00	up to -6
1.60	1.60 High Index	ADD 1.50 to 3.50	-10.00 to +8.00	up to -6
1.60	1.60 Transitions®	ADD 1.50 to 3.50	-10.00 to +8.00	up to -6
1.60	1.60 epik™ BLUE	ADD 1.50 to 3.50	-10.00 to +8.00	up to -6
1.67	1.67 High Index	ADD 1.50 to 3.50	-10.00 to +8.00	up to -6
1.67	1.67 Transitions®	ADD 1.50 to 3.50	-10.00 to +8.00	up to -6
1.67	1.67 epik™ BLUE	ADD 1.50 to 3.50	-10.00 to +8.00	up to -6

Fitting heights variable 14-35mm.

This lens design is optimized for patient Rx, is Rx Compensated, and has POW capabilities. Power ranges vary due to blank thickness and manufacturer. May also be affected by patient frame and Rx (prism).

epik™ iBifocal™

A breakthrough for D28 wearers - no line.

Index	Material Name	ADD	Power Ranges	CYL
1.50	CR39 Plastic	ADD .75 to 3.50	-9.25 to +7.50	up to -6
1.50	CR39 Transitions®	ADD .75 to 3.50	-9.50 to +7.50	up to -6
1.50	CR39 Polarized	ADD .75 to 3.50	-6.50 to +5.50	up to -6
1.50	CR39 epik™ BLUE	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.53	Trivex®	ADD .75 to 3.50	-9.00 to +6.00	up to -6
1.53	Trivex® Transitions®	ADD .75 to 3.50	-7.50 to +5.00	up to -4
1.53	Trivex® Polarized	ADD .75 to 3.50	-4.50 to +6.00	up to -4
1.59	PolyTrue™	ADD .75 to 3.50	-9.00 to +7.00	up to -6
1.59	Poly Transitions®	ADD .75 to 3.50	-9.50 to +7.50	up to -6
1.59	Poly Polarized	ADD .75 to 3.50	-10.00 to +7.00	up to -6
1.59	Poly epik™ BLUE	ADD .75 to 35.0	-9.00 to +7.00	up to -6
1.60	1.60 High Index	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.60	1.60 Transitions®	ADD .75 to 3.50	-10.00 to +8.00	up to -6
1.60	1.60 Polarized	ADD .75 to 3.50	-10.00 to +7.00	up to -6
1.60	1.60 epik™ BLUE	ADD .75 to 3.50	-10.00 to +8.00	up to -6

*Power ranges vary due to blank thickness and manufacturer.
May also be affected by patient frame and Rx (prism).*

ND4™ & epik™ DRV

Digital Device • Relaxed vision for pre-presbyopic patients.

Index	Material Name	Power Ranges	CYL
1.50	CR39 Plastic	-9.25 to +7.50	up to -6
1.50	CR39 Transitions®	-9.50 to +7.50	up to -6
1.50	CR39 Polarized	-6.50 to +5.50	up to -6
1.50	CR39 epik™ BLUE	-10.00 to +8.00	up to -6
1.53	Trivex®	-9.00 to +6.00	up to -6
1.53	Trivex® Transitions®	-7.50 to +5.00	up to -4
1.53	Trivex® Polarized	-4.50 to +6.00	up to -4
1.53	Trivex® epik™ BLUE	-8.00 to +6.00	up to -5
1.59	PolyTrue™	-9.00 to +7.00	up to -6
1.59	Poly Transitions®	-9.50 to +7.50	up to -6
1.59	Poly Polarized	-9.50 to +7.50	up to -6
1.59	Poly epik™ BLUE	-9.00 to +7.00	up to -6
1.60	1.60 High Index	-10.00 to +8.00	up to -6
1.60	1.60 Transitions®	-10.00 to +8.00	up to -6
1.60	1.60 Polarized	-10.00 to +7.00	up to -6
1.60	1.60 epik™ BLUE	-10.00 to +8.00	up to -6
1.67	1.67 High Index	-10.00 to +8.00	up to -6
1.67	1.67 Transitions®	-10.00 to +8.00	up to -6
1.67	1.67 Polarized	-10.00 to +8.00	up to -6
1.67	1.67 epik™ BLUE	-10.00 to +8.00	up to -6
1.74	1.74 High Index	-10.00 to +7.50	up to -6
1.74	1.74 Transitions®	-10.00 to +7.50	up to -6
1.74	1.74 epik™ BLUE	-10.00 to +7.50	up to -6

*Fit like a progressive, fitting heights variable 13-35mm.**This lens design is optimized for patient Rx.**Power ranges vary due to blank thickness and manufacturer.**May also be affected by patient frame and Rx (prism).***ND4™ & epik™ SV-Ultimate**

Ultra High Definition • Use to upgrade to all single vision.

Index	Material Name	Power Range	CYL
1.50	CR39 Plastic	-10.00 to +8.00	up to -6
1.50	CR39 Transitions®	-9.00 to +5.50	up to -6
1.50	CR39 Polarized	-8.00 to +7.00	up to -6
1.50	CR39 epik™ BLUE	-10.00 to +8.00	up to -6
1.53	Trivex®	-10.00 to +7.00	up to -6
1.53	Trivex® Transitions®	-7.50 to +6.50	up to -4
1.53	Trivex® Polarized	-7.00 to +7.50	up to -6
1.53	Trivex® epik™ BLUE	-8.00 to +6.00	up to -5
1.59	PolyTrue™	-10.00 to +8.00	up to -6
1.59	Poly Transitions®	-10.50 to +8.00	up to -6
1.59	Poly Polarized	-10.75 to +8.00	up to -6
1.59	Poly epik™ BLUE	-10.00 to +8.00	up to -6
1.60	1.60 High Index	-14.00 to +8.00	up to -6
1.60	1.60 Transitions®	-12.50 to +8.50	up to -6
1.60	1.60 Polarized	-11.50 to +8.00	up to -6
1.60	1.60 epik™ BLUE	-11.00 to +8.00	up to -6
1.67	1.67 High Index	-13.75 to +8.00	up to -6
1.67	1.67 Transitions®	-12.00 to +8.00	up to -6
1.67	1.67 Polarized	-12.00 to +9.00	up to -6
1.67	1.67 epik™ BLUE	-11.50 to +8.00	up to -6
1.74	1.74 High Index	-18.00 to +8.00	up to -6
1.74	1.74 Transitions®	-17.50 to +8.00	up to -6
1.74	1.74 epik™ BLUE	-18.00 to +8.00	up to -6

*For best performance specify Optical Center Height.**This lens design is optimized for patient RX, and has NO POW capabilities. Power ranges vary due to blank thickness and manufacturer. May also be affected by patient frame and Rx (prism).*

Stock Lens Ranges • CR39

CR-39 Factory Hard Coat

		CYL							
Power	SPH	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00
3.00									
2.75									
2.50									
2.25									
2.00									
1.75									
1.50									
1.25									
1.00									
0.75									
0.50									
+0.25									
0.00									
CYL →		-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00
0.00									
-0.25									
-0.50									
-0.75									
-1.00									
-1.25									
-1.50									
-1.75									
-2.00									
-2.25									
-2.50									
-2.75									
-3.00									
-3.25									
-3.50									
-3.75									
-4.00									
Power	SPH	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00

CR-39 with  NON-GLARE

		CYL							
Power	SPH	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00
3.00									
2.75									
2.50									
2.25									
2.00									
1.75									
1.50									
1.25									
1.00									
0.75									
0.50									
+0.25									
0.00									
CYL →		-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00
0.00									
-0.25									
-0.50									
-0.75									
-1.00									
-1.25									
-1.50									
-1.75									
-2.00									
-2.25									
-2.50									
-2.75									
-3.00									
-3.25									
-3.50									
-3.75									
-4.00									
Power	SPH	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00

Polarized Plano

includes factory hard coat
available in CR39: gray, brown, or G15

75 mm Blank Size

Polarized Plano Mirrors

includes factory hard coat
available in CR39: gray with blue, gold, or silver mirror

75 mm Blank Size

PLANO: Must order as Plano, do not select Single Vision under lens style when ordering. Only eligible if both lenses are Plano.

Stock Lens Ranges • PolyTrue™

PolyTrue™ Hard Coat

Power	SPH	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00	-2.25	-2.50	-2.75	-3.00	-3.25	-3.50	-3.75	-4.00
3.00																	
2.75																	
2.50																	
2.25																	
2.00																	
1.75																	
1.50																	
1.25																	
1.00																	
0.75																	
0.50																	
+2.5																	
0.00																	
CYL →		-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00	-2.25	-2.50	-2.75	-3.00	-3.25	-3.50	-3.75	-4.00
0.00																	
-0.25																	
-0.50																	
-0.75																	
-1.00																	
-1.25																	
-1.50																	
-1.75																	
-2.00																	
-2.25																	
-2.50																	
-2.75																	
-3.00																	
-3.25																	
-3.50																	
-3.75																	
-4.00																	
-4.25																	
-4.50																	
-4.75																	
-5.00																	
-5.25																	
-5.50																	
-5.75																	
-6.00																	
-6.25																	
-6.50																	
-6.75																	
-7.00																	
-7.25																	
-7.50																	
-7.75																	
-8.00																	
Power	SPH	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00	-2.25	-2.50	-2.75	-3.00	-3.25	-3.50	-3.75	-4.00

PolyTrue™ with *Infinity*
NON-GLARE

Power	SPH	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00	-2.25	-2.50	-2.75	-3.00	-3.25	-3.50	-3.75	-4.00
3.00																	
2.75																	
2.50																	
2.25																	
2.00																	
1.75																	
1.50																	
1.25																	
1.00																	
0.75																	
0.50																	
+2.5																	
0.00																	
CYL →		-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00	-2.25	-2.50	-2.75	-3.00	-3.25	-3.50	-3.75	-4.00
0.00																	
-0.25																	
-0.50																	
-0.75																	
-1.00																	
-1.25																	
-1.50																	
-1.75																	
-2.00																	
-2.25																	
-2.50																	
-2.75																	
-3.00																	
-3.25																	
-3.50																	
-3.75																	
-4.00																	
-4.25																	
-4.50																	
-4.75																	
-5.00																	
-5.25																	
-5.50																	
-5.75																	
-6.00																	
-6.25																	
-6.50																	
-6.75																	
-7.00																	
-7.25																	
-7.50																	
-7.75																	
-8.00																	
Power	SPH	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00	-2.25	-2.50	-2.75	-3.00	-3.25	-3.50	-3.75	-4.00

Stock Lens Ranges • PolyTrue™

PolyTrue™ Transitions® 8 Grey with **Infinity** NON-GLARE

		CYL							
Power	SPH	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00
0.00									
-0.25									
-0.50									
-0.75									
-1.00									
-1.25									
-1.50									
-1.75									
-2.00									
-2.25									
-2.50									
-2.75									
-3.00									
-3.25									
-3.50									
-3.75									
-4.00									
-4.25									
-4.50									
-4.75									
-5.00									
-5.25									
-5.50									
-5.75									
-6.00									
Power	SPH	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00

70 mm Blank Size

PolyTrue™ epik™ BLUE with **Infinity** NON-GLARE

		CYL							
Power	SPH	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00
4.00									
3.75									
3.50									
3.25									
3.00									
2.75									
2.50									
2.25									
2.00									
1.75									
1.50									
1.25									
1.00									
0.75									
0.50									
+0.25									
0.00									
CYL →		-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00
0.00									
-0.25									
-0.50									
-0.75									
-1.00									
-1.25									
-1.50									
-1.75									
-2.00									
-2.25									
-2.50									
-2.75									
-3.00									
-3.25									
-3.50									
-3.75									
-4.00									
-4.25									
-4.50									
-4.75									
-5.00									
-5.25									
-5.50									
-5.75									
-6.00									
-6.25									
-6.50									
-6.75									
-7.00									
-7.25									
-7.50									
-7.75									
-8.00									
Power	SPH	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00

63 mm Blank Size

76 mm Blank Size

epik™ BLUE Plano
includes Infinity® Non-Glare
available in poly and 1.60

75 mm Blank Size

PLANO: Must order as Plano, do not select Single Vision under lens style when ordering. Only eligible if both lenses are Plano.

Stock Lens Ranges • TRIVEX®

TRIVEX with *Infinity*
NONGLARE

		CYL							
Power	SPH	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00
4.00									
3.75									
3.50									
3.25									
3.00									
2.75									
2.50									
2.25									
2.00									
1.75									
1.50									
1.25									
1.00									
0.75									
0.50									
+0.25									
0.00									
CYL →		-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00
0.00									
-0.25									
-0.50									
-0.75									
-1.00									
-1.25									
-1.50									
-1.75									
-2.00									
-2.25									
-2.50									
-2.75									
-3.00									
-3.25									
-3.50									
-3.75									
-4.00									
-4.25									
-4.50									
-4.75									
-5.00									
-5.25									
-5.50									
-5.75									
-6.00									
-6.25									
-6.50									
-6.75									
-7.00									
Power	SPH	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00

65 mm Blank Size

70 mm Blank Size

TRIVEX Factory Hard Coat

		CYL							
Power	SPH	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00
4.00									
3.75									
3.50									
3.25									
3.00									
2.75									
2.50									
2.25									
2.00									
1.75									
1.50									
1.25									
1.00									
0.75									
0.50									
+0.25									
0.00									
CYL →		-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00
0.00									
-0.25									
-0.50									
-0.75									
-1.00									
-1.25									
-1.50									
-1.75									
-2.00									
-2.25									
-2.50									
-2.75									
-3.00									
-3.25									
-3.50									
-3.75									
-4.00									
-4.25									
-4.50									
-4.75									
-5.00									
-5.25									
-5.50									
-5.75									
-6.00									
-6.25									
-6.50									
-6.75									
-7.00									
Power	SPH	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00

65 mm Blank Size

70 mm Blank Size

Stock Lens Ranges • High Index

1.60 HI with **Infinity**
NON-GLARE

Power	SPH	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00
4.00									
3.75									
3.50									
3.25									
3.00									
2.75									
2.50									
2.25									
2.00									
1.75									
1.50									
1.25									
1.00									
0.75									
0.50									
0.25									
0.00									
CYL →		-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00

65 mm Blank Size

1.67 HI with **Infinity**
NON-GLARE

Power	SPH	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00
-1.00									
-1.25									
-1.50									
-1.75									
-2.00									
-2.25									
-2.50									
-2.75									
-3.00									
-3.25									
-3.50									
-3.75									
-4.00									
-4.25									
-4.50									
-4.75									
-5.00									
-5.25									
-5.50									
-5.75									
-6.00									
-6.25									
-6.50									
-6.75									
-7.00									
-7.25									
-7.50									
-7.75									
-8.00									
-8.25									
-8.50									
-8.75									
-9.00									
-9.25									
-9.50									
-9.75									
-10.00									
Power	SPH	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00

75 mm Blank Size

1.74 HI with **Infinity**
NON-GLARE

Power	SPH	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50
-3.00							
-3.25							
-3.50							
-3.75							
-4.00							
-4.25							
-4.50							
-4.75							
-5.00							
-5.25							
-5.50							
-5.75							
-6.00							
-6.25							
-6.50							
-6.75							
-7.00							
-7.25							
-7.50							
-7.75							
-8.00							
-8.25							
-8.50							
-8.75							
-9.00							
-9.25							
-9.50							
-9.75							
-10.00							
Power	SPH	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50

75 mm Blank Size
-6.25 to -10.00 70mm Blank

epik™ BLUE Plano
includes Infinity® Non-Glare
available in poly and 1.60

75 mm Blank Size

PLANO: Must order as Plano, do not select Single Vision under lens style when ordering.
Only eligible if both lenses are Plano.